

Short Lower Assembly, 8 - 12"

A1a	Left
A1b	Right

Footplate and mounting hardware not included.

Please Note: Short Lower Assemblies can only be used on 5 inch footplates. They cannot be used on 7 inch footplates.

Ref#	Torque Spec.
B1	22.5 – 27.5 N • m (16.6 – 20.3 ft • lb)

B1 Mounting Hardware
The lower extension requires a quantity of 2 screws per side to mount. A quantity of 4 must be selected.

Medium Lower Assembly, 10 - 16"

C1a	Left
C1b	Right

Footplate and mounting hardware not included.

Shown on a 7" x 8" footplate, but Medium Lowers can be used on any Version 4 footplates.

Lower Mounting Hardware

Ref#	Torque Spec.
E1	20 – 25 N • m (14.7 – 18.4 ft • lb)

Please Note: The Mounting Hardware is **NOT** included in the Lower Assembly and must be ordered individually.

B1 Footplate Mounting Hardware
The lower extension requires a quantity of 2 screws per side to mount. A quantity of 4 must be selected.

Long Lower Assembly, 13 - 19"

D1a	Left
D1b	Right

Footplate and mounting hardware not included.

Shown on a 7" x 11" footplate, but Long Lowers can be used on any Version 4 footplates.

Ref#	Torque Spec.
B1	22.5 – 27.5 N • m (16.6 – 20.3 ft • lb)