

Moldable Upper Support

See Owner's Manual for details



Padded Foot Support Sizing

All sizes are **inches** unless specified

* Dual always .5 narrower than described to fit between the leg rest frame.



	Width A	Depth B	Standard	Tall	For Dual with divider only	
			Height C		Divider Width D	Leg Trough E
SINGLE	7.5	10	8.5	12.5	-	-
12	11.5	10.5	8.5	12.5	-	-
14	13.5	10.5	8.5	12.5	-	-
16	15.5	10.5	8.5	12.5	2.5	5
18	17.5	10.5	8.5	12.5	2.5	6
20	19.5	10.5	8.5	12.5	2.5	7
22	21.5	10.5	8.5	12.5	2.5	8
24	23.5	10.5	8.5	12.5	2.5	9

Dual* (divided or not divided)