


(A1a)	Small 20"W x 20"D
(A1b)	Medium 22"W x 22"D
(A1c)	Large 24"W x 24"D



Complete Assembly


(B1)

To Joystick Receiver Mount*

Web Links



(WL1)

Joystick Receiver Mounts



Right Side Joystick Cutout Mounting

Left Side Joystick Cutout Mounting



Sizing Width Chart S = Small, M = Medium, L = Large	Arm to Arm Width											
	10"W	11"W	12"W	13"W	14"W	16"W	17"W	18"W	19"W	20"W	21"W	22"W
TRU Balance® Flip Back Arms	S	S	S	S	S	S	M	ML	ML	ML	L	L
Cane Mounted Flip Back	S	S	S	S	S	SM	M	ML	L	L	L	n/a
Flex T-Arms	S	S	S	S	S	SM	M	ML	L	L	L	n/a
Dual Post Flip Back	S	S	S	S	S	SM	M	ML	L	L	L	n/a
Single Post H/A Removable**	S	S	S	S	S	SM	SM	ML	ML	L	L	L