




# Pelvic Belts & Pads


- Help maintain a neutral pelvic position, which directly affects seated body alignment and distributes pressure appropriately.
- Increase postural stability, decreasing the effort needed to sit upright and facilitating upper extremity movement.
- Contribute to proper head positioning and neck alignment.
- Reduce negative postural tendencies that can lead to orthopedic asymmetries and other secondary medical complications.

## 2 Point Pad Sizes




	Part Number	A	B	C
XS	IFPB101	4.25	1.75	1.25
S	IFPB102	5.5	1.75	1.25
M	IFPB103	7.5	2.5	1.75
L	IFPB104	9.5	2.75	2.25
XL	IFPB105	11.5	2.75	2.25

## Belt Sizes



	A	B	C	Buckle	Pad Sizes
Single Pull	32	27	1	PB	NP / XS
Single Pull	32	27	1	SR	NP
Dual Pull	24	24	1	PB	XS / S
Single Pull	32	27	1.5	PB / SR	NP / M / L
Dual Pull	24	24	1.5	PB / SR	M / L
Single Pull	32	27	2	PB	NP / L / XL
Single Pull	32	27	2	SR	L / XL
Dual Pull	24	24	2	PB / SR	L / XL



A	B	C	D	Buckle	Sizes
32	27	20	1	PB	XS
32	27	20	1.5	PB	S / M / L

**PB:** Push Button **SR:** Side Release **NP:** Non-Padded