



# SP

- ▶ **New SIM G base, fluid resistant overlay top and sides**
- ▶ **TPU Heat-Sealed**
- ▶ **Low Profile**
- ▶ **Lightweight Cooling Gel Topper**




# GLACIAL

**Coccyx and Sacral Relief**

**Contoured to increase pelvic stability and contact area**

**Mildly contoured, high density molded foam**



**Gel Topper**

- Contoured base provides pelvic stability, lower extremity support, and increases overall surface contact area.
- Pelvic well provides a pre-ischial shelf to prevent posterior pelvic tilt. A silicone gel overlay provides additional skin protection and reduces shearing forces.
- Coolcore is the standard cushion cover material that has proven temperature and moisture control.
- Fire resistance: ISO 8191-1, ISO 8191-2, CAL 117-2013
- Biological evaluation/cytotoxicity ISO 10993-5, 10993-10
- Sub-trochanter support to increase lateral stability
- Maintains alignment of the femurs medially and laterally

**WEIGHT CAPACITY**

**450 lbs.**

**OVERALL WEIGHT**

**3½ - 5 lbs.**

**DIMENSIONS**

**from 10in (25.40cm) x 10in (25.40cm)  
to 22in (55.88cm) x 24in (60.96cm)**

SOF<sup>1</sup> available

**PROFILE**

**3½ in (8.89cm)**

**MSRPS**

**\$440 - \$595**

**HCPCS**

**E2603 | E2604**

1. Special Order Form
2. CoolCore® is standard on all cushions except the Simplicity G.
3. Special order sizes available for all cushions except the Glacial SP. Contact sales for additional assistance
4. Allowances listed are estimations for cushions billed in conjunction with Group 3 power chair bases using the KU modifier and may vary by state