

- ▶ **Pre-Ischial Shelf and Lateral Adductor build-up**
- ▶ **Mild Contoured Base with High Density Closed Cell posterior lateral foam inserts**
- ▶ **Extra Soft SunMate FRG Topper**



PREMIERE

WEIGHT CAPACITY

500 lbs.

OVERALL WEIGHT

3½ lbs.

DIMENSIONS

from 10in (25.40cm) x 10in (25.40cm)
to 22in (55.88cm) x 24in (60.96cm)

SOF¹ available

PROFILE

3½ in (7.62cm)

MSRPS

\$325 - \$475

HCPCS

E2605 | E2606

Coccyx and Sacral Relief

Contoured to increase pelvic stability and contact area



- Contoured base provides lower extremity support and increases overall surface contact area
- Additional high density foam inserts at the posterior lateral aspect of the cushion to promote offloading for increased pelvic positioning
- Posterior cutout provides relief at the coccyx and lower sacral region
- Top layer of Extra Soft SunMate FRG foam provides optimal protection through immersive properties and improves micro-climate
- CoolCore is the standard cushion cover material that has proven temperature and moisture control
- Fire resistance: ISO 8191-1, ISO 8191-2, CAL 117-2013
- Biological evaluation/cytotoxicity ISO 10993-5, 10993-10
- Sub-trochanter support to increase lateral stability
- Maintains alignment of the femurs medially and laterally

High density closed cell lateral foam inserts

1. Special Order Form
 2. CoolCore® is standard on all cushions except the Simplicity G.
 3. Special order sizes available for all cushions except the Glacial SP. Contact sales for additional assistance
 4. Allowances listed are estimations for cushions billed in conjunction with Group 3 power chair bases using the KU modifier and may vary by state