

- ▶ **New SIM G base**
- ▶ **SunMate FRG Extra Extra Soft**
- ▶ **Blue Topper, Low Profile**
- ▶ **Lightweight**

**SP**



# ZEN

**WEIGHT CAPACITY**

**400 lbs.**

**OVERALL WEIGHT**

**3½ lbs.**

**DIMENSIONS**

from 10in (25.40cm) x 10in (25.40cm)  
to 24in (60.96cm) x 24in (60.96cm)

SOF<sup>1</sup> available

**PROFILE**

**3in (7.62cm)**

**MSRPS**

**\$315 - \$425**

**HCPCS**

**E2603 | E2604**

**Coccyx and Sacral Relief**

**Contoured to increase pelvic stability and contact area**

**Mildly contoured, high density molded foam**



- Contoured base provides pelvic stability, lower extremity support, and increases overall surface contact area
- Posterior cutout provides relief at the coccyx and lower sacral region.
- Top layer of Extra Extra Soft SunMate FRG foam provides optimal skin protection through immersive properties and improves micro-climate
- CoolCore is the standard cushion cover material that has proven temperature and moisture control
- Fire resistance: ISO 8191-1, ISO 8191-2, CAL 117-2013
- Biological evaluation /cytotoxicity ISO 10993-5, 10993-10
- Sub-trochanter support to increase lateral stability
- Maintains alignment of the femurs medially and laterally

1. Special Order Form  
 2. CoolCore® is standard on all cushions except the Simplicity G.  
 3. Special order sizes available for all cushions except the Glacial SP. Contact sales for additional assistance  
 4. Allowances listed are estimations for cushions billed in conjunction with Group 3 power chair bases using the KU modifier and may vary by state