

ALL

POSITIONING



HEAD SOLUTIONS AP OVERVIEW

ALL POSITIONING MALLEABLE HEAD PAD

The **All Positioning (AP)** head support combines occipital and suboccipital support in one head pad. The AP offers a soft padded support for the occipital to rest on during tilt or recline and adjustable arms that can be bent to help provide more lateral suboccipital support. The AP discourages lateral flexion by improving head positioning. Supports swing away system if additional lateral support is required.

DUAL LAYERED FOAM MALLEABLE FRAME



WASHABLE & REMOVABLE

Standard stretchable
Neoprene cover



BASIC POSITIONING


Support to the **occipital** and **suboccipital** with same pad

OCCIPITAL SUPPORT


To **support** during **tilt/recline**

SUBOCCIPITAL SUPPORT

For **positioning** and **comfort**



A **malleable** frame that helps to adjust the pad contour to fit client's head contour for better comfort and support




AP OPTIONS

AP PAD SIZES



Cover can be **customized** with our variety of options at www.StealthProducts.com/upholstery




	I.N.	Description	Head Pad Width	Head Pad Height	Suboccipital Concave
A	AP150M	Metric Extra Small AP	7in (17.78cm)	6in (15.24cm)	1/2in (12.70mm)
	AP150	Extra Small AP	7in (17.78cm)	6in (15.24cm)	1/2in (12.70mm)
B	AP250M	Metric Small AP	8in (20.32cm)	6¼in (15.88cm)	3/4in (19.05mm)
	AP250	Small AP	8in (20.32cm)	6¼in (15.88cm)	3/4in (19.05mm)
C	AP350M	Metric Medium AP	9in (22.86cm)	6½in (16.51cm)	1in (2.54cm)
	AP350	Medium AP	9in (22.86cm)	6½in (16.51cm)	1in (2.54cm)
D	AP450M	Metric Large AP	10in (25.40cm)	7½in (19.05cm)	1¼in (3.18cm)
	AP450	Large AP	10in (25.40cm)	7½in (19.05cm)	1¼in (3.18cm)
E	AP550M	Metric Extra Large AP	13in (33.02cm)	8in (20.32cm)	1½in (3.81cm)
	AP550	Extra Large AP	13in (33.02cm)	8in (20.32cm)	1½in (3.81cm)

Only metric version is compatible with the new UniLink™ head solutions hardware

AP ADD-ONS SWING AWAY

EXTENDS **HEAD PAD**
POSITIONING & **MORE**



MORE **DETAILS** ABOUT
THE **SWING AWAY RECEIVER**
<https://stlpro.site/stealth-swingaway>



SAR2B SERIES
with Arm Height Adjustment

The metric AP is compatible with the **SAR2B** (Low Profile) swing away series while the original is compatible with the **SUS9** series. The receiver will allow a caretaker or technician to position any mini-joystick, switches or temporal pads where they are needed most.

The advantage of the SAR2B series over the SUS9 series is the arm height adjustment feature which has been introduced with the new lighter design.

Either to enhance positioning or to complement alternative control options, the Swing Away will let a caretaker or technician configure the Combo to offer solutions to any challenge.



SUS9 SERIES



ALTERNATIVE TECHNOLOGY

The **All Positioning** can be used as a support for tilt/recline. It can supplement with embedded switches on a tray driving control configuration, or any other configuration. It is compatible with our **i-Drive**® Advanced Alternative Control System and **i-Connect**™ sensor switches power router.

<https://stlpro.site/i-drive>

<https://stlpro.site/stealth-icconnect>