

# OCCIPITAL STANDARD PAD SELECTION

## STANDARD OCCIPITAL PADS




	I.N.	Description	Format	Total Width	Total Height	Contour Depth
A	530M	4x2 Occipital Pad	Curved	4in (10.16cm)	2½in (6.35cm)	¼in (6.35mm)
B	520M	4x3 Occipital Pad	Curved	4in (10.16cm)	3¼in (8.26cm)	¼in (6.35mm)
C	510M	5x2 Flat Occipital Pad	Flat	5in (12.70cm)	2in (5.08cm)	-
D	500M	7x4 Occipital Pad	Curved	7¼in (18.42cm)	4¼in (10.80cm)	1in (2.54cm)

## OTHER OCCIPITAL COMPATIBLE PADS



### CP06CCFPM Pediatric Comfort Plus

**Total Dimensions**  
6 ¼in (15.88cm) Wide  
4 ¼in (10.80cm) Height



### CB305M LARGE COMBO

**Middle Pad Dimensions**  
5in (12.70cm) Wide  
3¼in (8.26cm) Height

### CB315M SMALL COMBO

**Middle Pad Dimensions**  
5¾in (14.61cm) Wide  
4in (10.16cm) Height