

Suggested
Angle
Range

12°

Pad length
should **not**
go beyond
temple area

a Head Diameter
4 3/8in to 7 1/8in
11.11cm to 18.10cm


b Head Array
Between Pads
Diameter
4 1/2in to 7 1/4in
11.43cm to 18.41cm

c Space Between
Head and Pads
1/8in to 1/4in
0.32cm to 0.64cm




Curved Head Array

Approximate* Dimensions



Curved Design for a more ergonomic fit

400
450



	(A)	(a)	(B)	(b)	(e)	(f)
400	5 1/2in (13.97 cm)	6 1/2in (16.51 cm)	4 3/8in (11.11 cm)	5 1/2in (13.97 cm)	4 1/2in (11.83 cm)	2 1/4in (5.71 cm)
450	5 1/4in (14.83 cm)	5 27/32in (14.83 cm)	4 5/8in (11.75 cm)	5 9/16in (14.13 cm)	3 1/2in (8.89 cm)	1 1/2in (3.81 cm)

(C) Pad Thickness **1in** (2.54 cm)

(d) Hardware Depth **1in** (2.54 cm)

* Final dimensions can be affected by cushion materials, foam density and upholstery seams.