


Complete Upper Assemblies



Standard	
A1	Pair
B1	Single (L/R)
Long	
C1	Pair
D1	Single (L/R)

Lower Extension



Upper	Lower	Length (cm)
Standard	Standard	23 - 39
Standard	Long	33 - 44
Long	Standard	23 - 39
Long	Long	33 - 53


The maximum length can be no greater than 6.3 cm less than the seat to floor height (8.9 cm on the Q6 Edge 3 Stretto w/ TB3 seat). **Example:** A seat to floor height of 44.3 cm allows a maximum legrest length of 38 cm (35.4 cm on the Q6 Edge 3 Stretto w/ TB3 seat).

Long
H1b

A quantity of 2 must be selected for the pair.


Standard
H1a

Heel Loops



K1a Right
K1b Left


Calf Pad & Mounting Hardware




(L1) is included with the complete upper assemblies. (F1) and (G1) must be selected separately.

Replacement Actuator

Only parts with callouts are serviceable. For parts without callouts the complete assembly must be selected.




Footplates



Footplates (L / R)	Seat
I1a	10.5 x 17.5 cm 35.5 cm W & greater
I1b	12.4 x 17.5 cm 40.6 cm W & greater
I1c	14.8 x 17.5 cm 45.7 cm W & greater

A quantity of 2 must be selected for the pair.

Hangers



Pair w/ Mounting Hardware
J1 - J2