

Gel Properties

Non-Toxic
if consumed

Maintain Consistency
on temperatures
under 200°F (93°C)

Does Not Bottom Out
when weight is applied

Passed Flammability Testing
per RESNA WC-1
and international requirements

No Harmful Chemicals
passes CAL Prop 65, REACH and RoHS

Safe for the Environment
is a green recyclable material

Glycerin Mineral Oil
little to no odor
compared to gels in
the market

Twill Fabric
thin fabric coat
protecting the gel
both sides

Twill Fabric
65/35 polyester/cotton
woven fabric treated
with **non-hazardous**
dyes and chemicals

1/2in
1.27cm)

